





SOLUTIONS FOR OPHTHALMOLOGY KITS POSITIONING

Easy Fit System Baskets

Kit is placed in:











General technical features:

BASKET and LID

- Completely made of stainless steel AISI304
- rounded edges with smooth surfaces and absence of sharp edges, to ensure maximum safety for the user during handling and avoid any possibility of breakage of soft packaging
- Basket perforation opening: 4x4 mm holes on the base to allow better water penetration during the washing phases and 2x2 mm on the perimeter
- Cover perforation opening: 4x4 mm perforations
- predisposition for fixing the lid placed on both sides of the basket
- · Lid with ratchet closure system, practical and safe
- Low weight thanks to special production technology
- electropolished finish to prevent high sterilization temperatures of altering the colour and facilitate drying process and improve corrosion resistance
- Suitable for steam and low temperature sterilization

REMOVABLE INSTRUMENT HOLDER

- The basket is used in combination with a series of silicone supports for the stable and safe fixing of surgical instruments during the decontamination, washing, packaging, sterilization and transport processes
- Instrument holder made of 100% non-toxic green silicone (latex free)
- Support bar made of stainless steel AISI304
- T20 torx fixing screw (No. 2 per stand) made of stainless steel
- removable support to allow maximum adjustment according to the length of the instrument to be positioned
- Suitable for steam and low temperature sterilization

FIXING PIN

- Made in stainless steel AISI304
- T20 torx fixing screw (No. 1 per pin) made of stainless steel
- to place the ring instruments or lock the mat
- Suitable for steam and low temperature sterilization

CIRCULAR PERFORATED BOX

- made in stainless steel AISI304
- dim. Ø75 mm with central hole Ø4 mm for fixing
- fastening system (T20 torx screw + nut) made of stainless steel
- Suitable for steam and low temperature sterilization

FLEXYMAT MULTI-PERFORATION MAT

- Material: 100% non-toxic platinum silicone elastomer in green colour (Latex free)
- Alternating rows of conical tips (every 2 mm):
 - o high conical drill bits: Ø2x15h mm (Ø base), space between one and the other §
 - olow conical tips: Ø1x4h mm (Ø base)
- Drainage and steam circulation holes: 4.5x4.5 mm
- rounded base feet: Ø2,5x2h mm (Ø base)
- Mat height: 19 mm (H 2 mm, excluding feet and conical tips)
- Suitable for steam and low temperature sterilization

REMOVABLE SOFT LIPS SUPPORT FOR INSTRUMENT STOPS

- Material: 100% non-toxic green silicone (latex free)
- Soft lips support bar and fixing nuts made of stainless steel AISI304
- positioned on the cover to lock the instrumentation (both for removable supports and for mat)
- Suitable for steam and low temperature sterilization













Solution A: Grid with Removable Supports





Basket base + internal tray

Basket base

Kit placement solutions:

STANDARD SOLUTIONS P	ag. 4
CHALAZION P	ag. 7
CATARACT P	ag. 9
DALK+PKP P	ag. 23
RETINAL DETACHMENT P	ag. 24
DMEK P	ag. 25
TEAR DUCTS P	ag. 27
DSAEK P	ag. 28
PHACOEMULSIFIER P	ag. 29
SPECIAL INSTRUMENTS P	ag. 30
GLAUCOMA P	ag. 31
PALPEBRA P	ag. 34
MAJOR EYELID P	ag. 35
PLUMBING P	ag. 36
PTERYGIUM P	ag. 37
STRABISMUS P	ag. 38
STRICTUROTOMIA P	ag. 40
CORNEAL SUTURE P	ag. 42
PALPEBRA SUTURE P	ag. 42
BASIC TRANSPLANTS P	ag. 43
VITRECTOMY	ag. 44
MISCELLANEOUS SOLUTIONS P	ag. 48
SOLUTION FOR WASHING IN A WASHER-PISINFECTOR	ag. 49



Customized solutions





STANDARD SOLUTIONS

Art.


see table pag. 5

Description:

STAINLESS STEEL PERORATED SHEET BASKET COMPLETE WITH LID with HANDLE and REMOVABLE SILICONE SUPPORTS



Art. 71MICRO



Art. 71MICRO



Art. 77/2MICRO



Art. 77/3MICRO



Art. 71MICRO



Art. 77MICRO



Art. 77/1MICRO



Art. 77/1MICRO





Art.	Dim.	Components	Type of supports Silicone	Number of spaces that can be equipped
73MICRO	380x240x45h mm	basket + lid + 4 removable supports	2 stars + 1 tree + 1 comb	36 (2 lines)
73MICRO-TL	380x240x45h mm	Basket + flat lid + 4 removable supports + 2 soft lips	2 stars + 1 tree + 1 comb	36 (2 lines)
72MICRO	300x240x50h mm	basket + lid + 5 removable supports	2 stars + 1 tree + 2 comb	36 (4 lines)
72MICRO-TL	300x240x50h mm	Basket + flat lid + 5 removable supports + 3 soft lips	2 stars + 1 tree + 2 comb	36 (4 lines)
71MICRO	240x240x45h mm	basket + lid + 4 removable supports	2 stars + 1 tree + 1 comb	36 (2 lines)
71MICRO-TL	240x240x45h mm	Basket + flat lid + 4 removable supports + 2 soft lips	2 stars + 1 tree + 1 comb	36 (2 lines)
77MICRO	270x170x35h mm	Basket + flat lid + 2 removable supports	2 stars	18 (1 line)
77/1MICRO	270x170x35h mm	Basket + flat lid + 2 removable supports	1 star + 1 comb	18 (1 line)
77CLP	250x155x35h mm	Basket + flat lid + 3 removable supports	2 stars + 1 comb	18 (2 lines)
77/2MICRO	200x170x35h mm	basket + lid + 2 removable supports	1 star + 1 comb	14 (1 line)
77/3MICRO	200x170x35h mm	basket + lid + 2 removable supports	2 stars	14 (1 line)
63JKIT-1	200x80x35h mm	basket + lid + 2 removable supports	2 stars	3 (1 line)
63JKIT-2	200x80x35h mm	basket + lid + 2 removable supports	2 stars	4 (1 line)

Series 73, 72 & 71 - Specific Features:

- Complete with standard lid with recessed handle
- TL series complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- Recommended packaging: 1/2 US container 300x300 or soft packing





Others - Specific Features:

- series dim. 270x170 mm and dim. 250x155 complete with flat lid with central grip handle
- series dim. 200x170 mm complete with standard cover with recessed handle
- 63J series without grab handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- recommended packaging: MINI container 285x185 or soft packaging









SILICONE SUPPORT GEOMETRY

Quick and easy fastening of all the different instrument geometries. All housings allow the instrument to be positioned and to facilitate cleaning, washing and sterilization of the most contaminated area

STAR





CHRISTMAS TREE





COMB










Art. 71KIT-12

Description:

STAINLESS STEEL BASKET L.F. 240x240x50h mm COMPLETE WITH LID with HANDLE, 3 REMOVABLE SILICONE SUPPORTS and 3 PINS



Basket base (standard lid)

- Dimensions: 240x240x50h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 3 removable supports
- 3 stainless steel pins for positioning ring instruments
- Recommended packaging: 1/2 US container 300x300 or soft packing










CHALAZION Art. 77/2MICRO

Description:

STAINLESS STEEL BASKET L.F. 200x170x35h mm COMPLETE WITH LID with HANDLE and 2 REMOVABLE SILICONE SUPPORTS



Basket base (standard lid)

- Dimensions: 200x170x35h mm
- Complete with standard lid with recessed handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- 2 removable supports
- recommended packaging: MINI container series 285x185 or soft packaging










CATARACT Art. 71KIT-6

Description:

STAINLESS STEEL BASKET L.F. 240x240x45h mm COMPLETE WITH LID with HANDLE, 7 REMOVABLE SILICONE SUPPORTS and 3 PINS



Basket base (standard lid)

- Dimensions: 240x240x45h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 7 removable supports
- 3 stainless steel pins for positioning ring instruments
- Recommended packaging: 1/2 US container 300x300 or soft packing










CATARACT Art. 71/4RE

Description:

STAINLESS STEEL BASKET L.F. 240x240x41h mm COMPLETE WITH LID with HANDLE and SILICONE SOFT LIPS, 3 REMOVABLE SILICONE SUPPORTS, MAT and MINI EQUIPPED INTERNAL BASKET



Basket base (standard lid)

- Dimensions: 240x240x41h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 3 removable supports
- 1 multi-perforated silicone mat
- 1 removable silicone soft lip in the lid
- 1 mini basket dim. 170x50x25h mm with hinged lid, mat and 2 silicone soft lips
- Recommended packaging: 1/2 US container 300x300 or soft packing









CATARACT Art. 71PSM

Description:

STAINLESS STEEL BASKET L.F. 240x240x70h mm COMPLETE WITH LID with HANDLE, 15 REMOVABLE SILICONE SUPPORTS, MAT and 2 PINS





Basket base + lid with supports for phacoemulsifier fixing

- Dimensions: 240x240x70h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 15 removable mounts
- 1 multi-perforated silicone mat
- 2 stainless steel pins to lock the mat
- Recommended packaging: 1/2 US container 300x300 or soft packing









CATARACT Art. 71KIT-21

Description:

STAINLESS STEEL BASKET L.F. 240x240x50h mm COMPLETE WITH LID with HANDLE, 4 REMOVABLE SILICONE SUPPORTS, 4 PINS and 2 MATS



Basket base (standard lid)

- Dimensions: 240x240x50h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 4 removable supports
- 4 stainless steel pins for positioning ring instruments
- 2 multi-perforated silicone mats
- Recommended packaging: 1/2 US container 300x300 or soft packing










CATARACT Art. 71PSM-1

Description:

STAINLESS STEEL BASKET L.F. 240x240x70h mm COMPLETE WITH LID with HANDLE and 7 REMOVABLE SILICONE SUPPORTS



Basket base (standard lid)

- Dimensions: 240x240x70h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- · Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 7 removable supports
- Recommended packaging: 1/2 US container 300x300 or soft packing









CATARACT Art. 71KIT-24

Description:

STAINLESS STEEL BASKET L.F. 240x240x50h mm COMPLETE WITH LID with HANDLE, 2 REMOVABLE SILICONE SUPPORTS, 3 PINS, SILICONE SOFT LIPS, CANNULA BOX and DIVIDER WITH STOPS



Basket base (standard lid)

- Dimensions: 240x240x50h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- · 2 removable supports
- 3 stainless steel pins for positioning ring instruments
- 1 removable silicone soft lip in the lid
- Circular perforated box for cannulas
- stainless steel partition L=200 mm complete with 2 stops
- Recommended packaging: 1/2 US container 300x300 or soft packing











CATARACT Art. 71KIT-20

Description:

STAINLESS STEEL BASKET L.F. 240x240x70h mm COMPLETE WITH LID with HANDLE, 4 REMOVABLE SILICONE SUPPORTS, MAT and 2 PINS





Basket base (standard lid)

- Dimensions: 240x240x70h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 4 removable supports
- 1 multi-perforated silicone mat
- 2 stainless steel pins to lock the mat
- Recommended packaging: 1/2 US container 300x300 or soft packing










CATARACT Art. 71KIT-27

Description:

STAINLESS STEEL BASKET L.F. 240x240x41h mm COMPLETE WITH LID with HANDLE, 2 REMOVABLE SILICONE SUPPORTS and SILICONE SOFT LIPS





Basket base

Cover with silicone soft lip

- Dimensions: 240x240x41h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 2 removable supports
- 1 removable silicone soft lip in the lid
- Recommended packaging: 1/2 US container 300x300 or soft packing










CATARACT Art. 71/1-CATARACT-SET

Description:

STAINLESS STEEL BASKET L.F. 240x240x50h mm COMPLETE WITH LID with HANDLE, 3 REMOVABLE SILICONE SUPPORTS and 12 PINS



Basket base (standard lid)

- Dimensions: 240x240x50h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 3 removable supports
- 16 stainless steel pins for positioning ring instruments and capsules
- Recommended packaging: 1/2 US container 300x300 or soft packing










CATARACT Art. 77KIT-1

Description:


STAINLESS STEEL BASKET L.F. 250x155x35h mm COMPLETE WITH LID with HANDLE, 5 REMOVABLE SILICONE SUPPORTS, MAT, SILICONE SOFT LIPS and 5 PINS



Basket base (standard lid with SILICONE SOFT LIP)

- Dimensions: 250x155x35h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- 5 removable supports
- 5 stainless steel pins for positioning the ring tools and locking the mat
- 1 multi-perforated silicone mat
- 1 removable silicone soft lip in the lid
- recommended packaging: MINI container series 285x185 or soft packaging










CATARACT BASIC TOOLS

Art. 77KIT-21

Description:

STAINLESS STEEL BASKET L.F. 270x170x35h mm COMPLETE WITH LID with HANDLE and 4 REMOVABLE SILICONE SUPPORTS



Basket base (standard lid)

- Dimensions: 270x170x35h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- 4 removable supports
- recommended packaging: MINI container series 285x185 or soft packaging









CATARACT + 2 PHATOEMULSIFIER

71KIT-1 Art.

Description:

STAINLESS STEEL BASKET L.F. 240x240x90h mm COMPLETE WITH LID with HANDLE and SILICONE SOFT LIP, INTERNAL TRAY, 11 REMOVABLE SILICONE SUPPORTS, MAT and 6 PINS



Basket base for 2 phacoemulsifiers



Internal tool tray



Base + Inner Tray



Lid with silicone soft lip

- Dimensions: 240x240x90h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid and without internal tray
- Upper inner tray made of stainless steel AISI304
- 11 removable supports
- 6 stainless steel pins for positioning the ring tools and locking the mat
- 1 multi-perforated silicone mat
- 1 removable silicone soft lip in the lid
- Recommended packaging: 1/2 US container 300x300 or soft packing










CATARACT SUPPLEMENT

Art.

71KIT-13

Description:

STAINLESS STEEL BASKET L.F. 240x240x45h mm COMPLETE WITH LID with HANDLE, 4 REMOVABLE SILICONE SUPPORTS and BOX FOR CANNULAS



Basket base (standard lid)

- Dimensions: 240x240x45h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 4 removable supports
- Circular perforated box for cannulas
- Recommended packaging: 1/2 US container 300x300 or soft packing










CATARACT SUPPLEMENT + PHATOEMULSIFIER


Art.

71KIT-23

Description:

STAINLESS STEEL BASKET L.F. 240x240x45h mm COMPLETE WITH LID with HANDLE, 10 REMOVABLE SILICONE SUPPORTS and BOX FOR CANNULAS





Basket base + lid with supports for phacoemulsifier fixing

- Dimensions: 240x240x45h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 10 removable supports
- · Circular perforated box for cannulas
- Recommended packaging: 1/2 US container 300x300 or soft packing









DALK+PKP Art. 71KIT-17

Description:

STAINLESS STEEL BASKET L.F. 240x240x45h mm COMPLETE WITH LID with HANDLE, 4 REMOVABLE SILICONE SUPPORTS, PIN and 2 BOXES FOR CANNULAS



Basket base (standard lid)

- Dimensions: 240x240x45h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- · Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 4 removable supports
- 2 circular perforated boxes for cannulas
- 1 stainless steel pin
- Recommended packaging: 1/2 US container 300x300 or soft packing










RETINAL DETACHMENT


Art.

71/3-RETINAL-SET

Description:

STAINLESS STEEL BASKET L.F. 240x240x90h mm COMPLETE WITH LID with HANDLE, INTERNAL TRAY, 3 REMOVABLE SILICONE SUPPORTS and 15 PINS





Basket base (standard lid)

Upper inner tray

- Dimensions: 240x240x90h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid and without internal tray
- Upper inner tray made of stainless steel AISI304
- 3 removable supports
- 15 stainless steel pins for positioning ring instruments and capsules
- Recommended packaging: 1/2 US container 300x300 or soft packing










DMEK Art. 71KIT-15

Description:

STAINLESS STEEL BASKET L.F. 240x240x70h mm COMPLETE WITH LID with HANDLE, 4 REMOVABLE SILICONE SUPPORTS, MINI BASKET and BOX FOR CANNULAS



Basket base (standard lid)

- Dimensions: 240x240x70h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 4 removable supports
- Mini basket with cannula holder support art. SCO4/M-H40 (marked: A-B-DX-SX)
- Circular perforated box for cannulas
- Recommended packaging: 1/2 US container 300x300 or soft packing











DMEK Art. 71KIT-22

Description:

STAINLESS STEEL BASKET L.F. 240x240x90h mm + STAINLESS STEEL BASKET L.F. 240x240x50h mm COMPLETE WITH LID with HANDLE, 7 REMOVABLE SILICONE SUPPORTS, MINI BASKET, CANNULA BOX and 3 PINS





Basket base 240x240x90, 2 removable silicone supports, mini-basket, cannula box and 3 pins

Basket base 240x240x50 complete with lid with handle and 5 removable silicone supports

- 1 basket dimensions: 240x240x90h mm with grab handles
- 1 basket dimensions: 240x240x50h mm without grab handles
- 1 flat lid with central grip handle
- Handle made of stainless steel AISI304
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 7 removable supports
- mini basket
- · Circular perforated box for cannulas
- 3 stainless steel pins for positioning ring instruments
- Recommended packaging: 1/2 US container 300x300 or soft packing











TEAR DUCTS Art. 75KIT-8

Description:

STAINLESS STEEL BASKET L.F. 540x240x90h mm COMPLETE WITH LID with HANDLE, INTERNAL TRAY, 3 DIVIDERS, MINI BASKET, 8 REMOVABLE SILICONE SUPPORTS, 6 PINS and CANNULAS BOX



Basket base (standard lid)



Upper inner tray

- Dimensions: 540x240x90h mm
- Complete with flat lid with central grip handle
- handles made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid and without internal tray
- Upper inner tray made of stainless steel AISI304
- 8 removable supports
- 6 stainless steel pins for positioning ring instruments
- 3 partitions positioned in the base of the basket
- mini basket
- Circular perforated box for cannulas
- Recommended packaging: 1 US 600x300 container or soft packing








DSAEK Art. 71KIT-16

Description:

STAINLESS STEEL BASKET L.F. 240x240x70h mm COMPLETE WITH LID with HANDLE, 4 REMOVABLE SILICONE SUPPORTS, MINI BASKET and BOX FOR CANNULAS



Basket base (standard lid)

- Dimensions: 240x240x70h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 4 removable supports
- Mini basket with cannula holder support art. SCO4-H40
- Circular perforated box for cannulas
- Recommended packaging: 1/2 US container 300x300 or soft packing









PHACOEMULSIFIER Art. 76E-PHACO-SET

Description:

STAINLESS STEEL BASKET L.F. 280x110x40h mm COMPLETE WITH LID, 4 REMOVABLE SILICONE SUPPORTS and 4 PINS



Basket base and lid

- Dimensions: 280x110x40h mm
- Complete with flat lid without grab handle
- Basket base without grab handles
- 4 removable supports
- 4 pins for cable positioning
- recommended packaging: MINI container series 300x130 or soft packaging










SPECIAL INSTRUMENTS


Art.

71KIT-18

Description:

STAINLESS STEEL BASKET L.F. 240x240x90h mm COMPLETE WITH LID with HANDLE, INTERNAL TRAY, MINI BASKET, 5 REMOVABLE SILICONE SUPPORTS, 3 PINS and BOX FOR **CANNULAS**






Basket base (standard lid)

Upper inner tray

- Dimensions: 240x240x90h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid and without internal tray
- Upper inner tray made of stainless steel AISI304
- 5 removable supports
- 3 stainless steel pins for positioning ring instruments
- Mini basket with cannula holder support art. SCO4-H40
- Circular perforated box for cannulas
- Recommended packaging: 1/2 US container 300x300 or soft packing











GLAUCOMA Art. 71KIT-8

Description:

STAINLESS STEEL BASKET L.F. 240x240x70h mm COMPLETE WITH LID with HANDLE, INTERNAL TRAY, 6 REMOVABLE SILICONE SUPPORTS, 3 PINS and BOX FOR CANNULAS





Basket base (standard lid)

Upper inner tray

- Dimensions: 240x240x70h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid and without internal tray
- Upper inner tray made of stainless steel AISI304
- 6 removable supports
- 3 stainless steel pins for positioning ring instruments
- Circular perforated box for cannulas
- Recommended packaging: 1/2 US container 300x300 or soft packing










GLAUCOMA Art. 77/1MICRO

Description:

STAINLESS STEEL BASKET L.F. 270x170x35h mm COMPLETE WITH LID with HANDLE and 2 REMOVABLE SILICONE SUPPORTS



Basket base (standard lid)

- Dimensions: 270x170x35h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- 2 removable supports
- recommended packaging: MINI container series 285x185 or soft packaging










GLAUCOMA BASE KIT

Art.

77KIT-21

Description:

STAINLESS STEEL BASKET L.F. 270x170x35h mm COMPLETE WITH LID with HANDLE and 3 REMOVABLE SILICONE SUPPORTS



Basket base (standard lid)

- Dimensions: 270x170x35h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- 3 removable supports
- recommended packaging: MINI container series 285x185 or soft packaging










EYELID Art. 71/1-LID-SET

Description:

STAINLESS STEEL BASKET L.F. 240x240x50h mm COMPLETE WITH LID with HANDLE, 2 REMOVABLE SILICONE SUPPORTS and 12 PINS



Basket base (standard lid)

- Dimensions: 240x240x50h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 2 removable supports
- 12 stainless steel pins for positioning the ring instruments of the capsules
- Recommended packaging: 1/2 US container 300x300 or soft packing










MAJOR EYELID Art. 71KIT-9

Description:

STAINLESS STEEL BASKET L.F. 240x240x50h mm COMPLETE WITH LID with HANDLE, 2 REMOVABLE SILICONE SUPPORTS and 3 PINS



Basket base (standard lid)

- Dimensions: 240x240x50h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- · stackable baskets without lid
- 2 removable supports
- 3 stainless steel pins for positioning ring instruments
- Recommended packaging: 1/2 US container 300x300 or soft packing











PLUMBING Art. 71KIT-7

Description:

STAINLESS STEEL BASKET L.F. 240x240x70h mm COMPLETE WITH LID with HANDLE, INTERNAL TRAY, MINI BASKET, 8 REMOVABLE SILICONE SUPPORTS and 3 PINS





Basket base (standard lid)

Upper inner tray

- Dimensions: 240x240x70h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- · Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid and without full tray
- Upper inner tray made of stainless steel AISI304
- 8 removable supports
- 3 stainless steel pins for positioning ring instruments
- mini basket
- Recommended packaging: 1/2 US container 300x300 or soft packing










PTERYGIUM Art. 77KIT-15

Description:


STAINLESS STEEL BASKET L.F. 270x170x35h mm COMPLETE WITH LID with HANDLE, 2 REMOVABLE SILICONE SUPPORTS and 3 PINS



Basket base (standard lid)

- Dimensions: 270x170x35h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- 2 removable supports
- 3 stainless steel pins for positioning ring instruments
- recommended packaging: MINI container series 285x185 or soft packaging











STRABISMUS Art. 71KIT-11

Description:

STAINLESS STEEL BASKET L.F. 240x240x70h mm COMPLETE WITH LID with HANDLE, INTERNAL TRAY, MINI BASKET, 2 REMOVABLE SILICONE SUPPORTS and 6 PINS





Basket base (standard lid)

Upper inner tray

- Dimensions: 240x240x70h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid and without internal tray
- Upper inner tray made of stainless steel AISI304
- 2 removable supports
- 6 stainless steel pins for positioning ring instruments
- mini basket
- Recommended packaging: 1/2 US container 300x300 or soft packing










STRABISMUS Art. 71/3-SQUINT-SET

Description:

STAINLESS STEEL BASKET L.F. 240x240x90h mm COMPLETE WITH LID with HANDLE, INTERNAL TRAY, 5 REMOVABLE SILICONE SUPPORTS and 16 PINS





Basket base (standard lid)

Upper inner tray

- Dimensions: 240x240x90h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid and without internal tray
- Upper inner tray made of stainless steel AISI304
- 5 removable supports
- 16 stainless steel pins for positioning ring instruments and capsules
- Recommended packaging: 1/2 US container 300x300 or soft packing










STRICTUROTOMIA Art. 77KIT-18

Description:


STAINLESS STEEL BASKET L.F. 270x170x35h mm COMPLETE WITH LID with HANDLE, 2 REMOVABLE SILICONE SUPPORTS, 3 PINS and BOX FOR CANNULAS



Basket base (standard lid)

- Dimensions: 270x170x35h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- 2 removable supports
- 3 stainless steel pins for positioning ring instruments
- Circular perforated box for cannulas
- recommended packaging: MINI container series 285x185 or soft packaging










STRICTUROTOMIA Art. 71KIT-19

Description:

STAINLESS STEEL BASKET L.F. 240x240x41h mm COMPLETE WITH LID with HANDLE, 4 REMOVABLE SILICONE SUPPORTS, 3 PINS and BOX FOR CANNULAS



Basket base (standard lid)

- Dimensions: 240x240x41h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 4 removable supports
- 3 stainless steel pins for positioning ring instruments
- · Circular perforated box for cannulas
- Recommended packaging: 1/2 US container 300x300 or soft packing










CORNEAL SUTURE - EYELID SUTURE

Art.

77KIT-16

Description:

STAINLESS STEEL BASKET L.F. 270x170x35h mm COMPLETE WITH LID with HANDLE, 2 REMOVABLE SILICONE SUPPORTS and 3 PINS



Basket base (standard lid)

- Dimensions: 270x170x35h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- 2 removable supports
- 3 stainless steel pins for positioning ring instruments
- recommended packaging: MINI container series 285x185 or soft packaging










TRANSPLANT BASE Art. 71KIT-14

Description:

STAINLESS STEEL BASKET L.F. 240x240x70h mm COMPLETE WITH LID with HANDLE, 4 REMOVABLE SILICONE SUPPORTS and 3 PINS



Basket base (standard lid)

- Dimensions: 240x240x70h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 4 removable supports
- 3 stainless steel pins for positioning ring instruments
- Recommended packaging: 1/2 US container 300x300 or soft packing











VITRECTOMY Art. 71KIT-10

Description:

STAINLESS STEEL BASKET L.F. 240x240x70h mm COMPLETE WITH LID with HANDLE, INTERNAL TRAY, 5 REMOVABLE SILICONE SUPPORTS, 3 PINS and BOX FOR CANNULAS





Basket base (standard lid)

Upper inner tray

- Dimensions: 240x240x70h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid and without internal tray
- Upper inner tray made of stainless steel AISI304
- 5 removable supports
- 3 stainless steel pins for positioning ring instruments
- Circular perforated box for cannulas
- Recommended packaging: 1/2 US container 300x300 or soft packing










VITRECTOMY Art. 71KIT-4

Description:

STAINLESS STEEL BASKET L.F. 240x240x41h mm COMPLETE WITH LID with HANDLE, 2 REMOVABLE SILICONE SUPPORTS, 3 PINS, BOX FOR CANNULAS and SILICONE SOFT LIP





Basket base

Cover with silicone soft lip

- Dimensions: 240x240x41h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- 2 removable supports
- 3 stainless steel pins for positioning ring instruments
- Circular perforated box for cannulas
- 1 removable silicone soft lip in the lid
- Recommended packaging: 1/2 US container 300x300 or soft packing











VITRECTOMY Art. 73KIT-9

Description:

STAINLESS STEEL BASKET L.F. 380x240x90h mm COMPLETE WITH LID with HANDLE, INTERNAL TRAY, 8 REMOVABLE SILICONE SUPPORTS and 3 PINS



Basket base (standard lid)



Upper inner tray

- Dimensions: 380x240x90h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid and without internal tray
- Upper inner tray made of stainless steel AISI304
- 8 removable supports
- 3 stainless steel pins for positioning ring instruments
- Recommended packaging: 3/4 US container 460x300 or soft packing









VITRECTOMY Art. 73KIT-10

Description:

STAINLESS STEEL BASKET L.F. 380x240x41h mm COMPLETE WITH LID with HANDLE, 5 REMOVABLE SILICONE SUPPORTS, MAT, 2 SILICONE SOFT LIPS, 5 PINS and CANNULAS BOX



Basket base (standard lid)

- Dimensions: 380x240x41h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid and without internal tray
- 5 removable supports
- 1 multi-perforated silicone mat
- 2 removable silicone soft lips in the lid (1xW=220 mm + 1xW=130 mm)
- 5 stainless steel pins for positioning the ring tools and locking the mat
- Circular perforated box for cannulas
- Recommended packaging: 3/4 US container 460x300 or soft packing









MISCELLANEOUS SOLUTIONS

Art.

See table below

Art.	Description	Instrument	lmage
67SMKIT-1	STAINLESS STEEL BASKET L.F. 170x135x35h mm COMPLETE WITH FLAT LID without HANDLE and 4 REMOVABLE SILICONE SUPPORTS	2 handles for coaxial cannulas	
63KIT-4	STAINLESS STEEL BASKET L.F. 170x50x30h mm COMPLETE WITH FLAT LID without HANDLE and 2 REMOVABLE SILICONE SUPPORTS	1 handle for coaxial cannulas	
63KIT-6	STAINLESS STEEL BASKET L.F. 170x50x30h mm COMPLETE WITH FLAT LID without HANDLE and 2 REMOVABLE SILICONE SUPPORTS	1 tool	
SCO4-H40	STAINLESS STEEL SUPPORT with SILICONE RETAINER OPHTHALMIC CANNULA HOLDER dim. 55x55x37h	4 CANNULAS h max 40 mm	
SCO4/M-H40	STAINLESS STEEL SUPPORT with SILICONE RETAINER OPHTHALMIC CANNULA HOLDER dim. 55x55x37h with A-B-DX-SX marking	4 CANNULAS h max 40 mm	A SI SX
KIT-S5	KIT 6 REMOVABLE SILICONE SUPPORTS FOR FIXING FACO ON THE LID (media position: cover)	1 PHACOEMULSIFIER	






WASHER-DISINFECTOR WASHING SOLUTION

Art.

77KIT-22


Description:

STAINLESS STEEL BASKET L.F. 270x170x35h mm COMPLETE WITH LID with HANDLE, 4 REMOVABLE SILICONE SUPPORTS, MAT, 2 SILICONE SOFT LIPS and 2 PINS





Basket base + lid with silicone soft lip



Removable silicone mounts to lock 2 cannulated instruments



Basket base equipped with PASSAGE HOLES FOR WASHER-DISINFECTOR WASHING TUBES (2 holes for each corner of the basket, total 8 holes)

- Dimensions: 270x170x35h mm
- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base without grab handles
- 4 removable supports
- 1 multi-perforated silicone mat
- 2 stainless steel pins to lock the mat
- 2 silicone soft lips
- Basket base with 8 passage holes for washer-disinfector washing pipes
- if used for sterilization, recommended packaging: MINI container series 285x185 or soft packaging











Solution B: Grill with Mat and silicone soft lip in the lid









SINGLE SHELF SOLUTIONS

Art.

See table below


Description:L.F. STAINLESS STEEL BASKET COMPLETE WITH LID with HANDLE, MULTI-PERFORATED SILICONE MAT and SILICONE SOFT LIPS IN THE LID

Art.	Dim.	Components	Features	Number of instruments *
75/4TL	540x240x41h mm	basket + lid + Mat + 5 silicone soft lips	Flat lid and removable silicone soft lips	60
74/4TL	480x240x41h mm	basket + lid + Mat + 4 silicone soft lips	Flat lid and removable silicone soft lips	48
74/4TLQ	480x240x41h mm	basket + lid + Mat + 2 silicone soft lips	flat lid and fixed silicone soft lips L=440 mm	48
73/4TL	380x240x41h mm	basket + lid + Mat + 2 silicone soft lips	Flat lid and removable silicone soft lips	36
73/4TLQ	380x240x41h mm	basket + lid + Mat + 2 silicone soft lips	Flat lid and fixed silicone soft lips	36
71/4TL	240x240x41h mm	basket + lid + Mat + 2 silicone soft lips	Flat lid and removable silicone soft lips	24
71/4TLQ	240x240x41h mm	basket + lid + Mat + 2 silicone soft lips	Flat lid and fixed silicone soft lips	24
76TL	280x110x40h mm	basket + lid + Mat + 2 silicone soft lips	lid with recessed handle and removable silicone soft lips	6





" Silicone soft lips " position:




2719T	270x190x35h mm	basket + lid + Mat + 2 longitudinal silicone soft lips	lid with recessed handle and fixed silicone soft lips	12
77TH	270x170x35h mm	basket + lid + Mat + 2 longitudinal silicone soft lips	lid with recessed handle and fixed silicone soft lips	12
77TH-P	270x170x35h mm	basket + lid + Mat + 2 longitudinal silicone soft lips	Flat lid with handle and fixed silicone soft lips	12
77/2TH	200x170x35h mm	basket + lid + Mat + 2 longitudinal silicone soft lips	lid with recessed handle and fixed silicone soft lips	8
63MDL	220x50x35h mm	basket + lid + Mat + 2 transverse silicone soft lips	handleless lid and fixed silicone soft lips	8
67M	170x135x35h mm	basket + lid + Mat + 2 transverse silicone soft lips	handleless lid and fixed silicone soft lips	6
63T	170x50x25h mm	basket + lid + Mat + 2 transverse silicone soft lips	handleless lid and fixed silicone soft lips	1

^{*} the number of tools that can be placed within the grid is indicative and changes according to their geometry









Set > 24 Irons - Specifications:

- Complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- Basket base with grab handles
- Multi-perforation silicone mat
- Removable silicone soft lips
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement, allows stacking and a smoother load on the washer-disinfector trolleys
- stackable baskets without lid
- Recommended packaging: 1/2 US, 3/4 US or 1 US container 300/460/600x300 (depending on the size
 of the grid) or soft packaging





Set ≤ 12 irons - Specific Features:


- Complete with standard lid with recessed handle
- art. 77TH-P complete with flat lid with central grip handle
- Handle made of stainless steel AISI304
- art. 63MDL, 67M and 63T without grab handle
- Basket base without grab handles
- Multi-perforation silicone mat
- silicone soft lips fixed in the lid
- art. 76TL with removable silicone soft lips
- recommended packaging: MINI container series 285x185 or soft packaging









Art. 77/2TH





Art. 74/4TL



Art. 73/4TL







DOUBLE & SINGLE SHELF SOLUTIONS

Art.

See table below

Description:

L.F. STAINLESS STEEL BASKET COMPLETE WITH HINGED LID with HANDLE, 2 MULTI-PERFORATED SILICONE MATS, INTERNAL SHELF and 4 SILICONE SOFT LIPS (2 LID and 2 INTERNAL SHELF)







Art. 71/2-WMS

lid NOT detachable from the base

Art.	Dim.	Components	Features	Number of instruments
71/2-WMS	DOUBLE SHELF 240x240x60h mm	basket + hinged lid + inner shelf + 2 Mats + 4 silicone soft lips	flat lid that cannot be removed from the base, removable shelf, removable silicone soft lips	48
71/5C-WMS	SINGLE SHELF 240x240x35h mm	basket + hinged lid + mat + 2 silicone soft lips	flat lid that cannot be detached from the base, removable silicone soft lips	24
77/2-WMS	DOUBLE SHELF 270x170x60h mm	basket + hinged lid + inner shelf + 2 Mats + 4 silicone soft lips	flat lid that cannot be removed from the base, removable shelf, removable silicone soft lips	24
77THC-WMS	SINGLE SHELF 270x170x35h mm	basket + hinged lid + mat + 2 silicone soft lips	flat lid that cannot be detached from the base, removable silicone soft lips	12

^{*} the number of tools that can be placed within the grid is indicative and varies according to their geometry









Art. 77/2-WMS

Features Specifications dim. 240x240:

- Complete with flat lid with central grip handle. The lid is hinged and cannot be detached from the base
- Handle made of stainless steel AISI304
- Basket base without grab handles
- stainless steel wire foot (h = 6 mm) positioned at the base of the basket. The foot acts as a reinforcement and allows a smoother load on the washer-disinfector trolleys
- "double shelf" version with removable internal shelf made of stainless steel AISI304
- Multi-perforated silicone mats
- removable silicone soft lips (in the "double shelf" version there are 2 on the lid and 2 on the bottom shelf)
- Recommended packaging: 1/2 US container 300x300 or soft packing





Features Specifications dim. 270x170:

- Complete with flat lid with central grip handle. The lid is hinged and cannot be detached from the base
- Handle made of stainless steel AISI304
- Basket base without grab handles
- "double shelf" version with removable internal shelf made of stainless steel AISI304
- Multi-perforated silicone mats
- removable silicone soft lips instrument (in the "double shelf" version there are 2 on the lid and 2 on the bottom shelf)
- recommended packaging: MINI container series 285x185 or soft packaging









Producer:	C.B.M. SRL
Product Origin:	Made in Italy
Conformity:	Regulation 2017/745/EU "Medical Devices", class I
Quality Standards:	UNI IN ISO 13485:2021
Guarantee:	24 months
Packaging Type:	container or soft packing
Sterilization mode:	Steam & Low Temperature
Pieces per pack:	1
Primary packaging:	Single polyethylene bag or single carton box with information label for packaging

Use:

The perforated sheet metal baskets complete with internal accessories are used for the positioning of ophthalmology instruments. They adapt perfectly to every need and are essential during the various phases of the sterilization process (decontamination, washing (ultrasonics and washer-disinfectors), packaging, sterilization, handling and storage of surgical devices. In fact:

- are entirely made from perforated stainless steel sheet AISI304
- Thanks to the smooth surfaces, the absence of edges and the rounded corners, the basket guarantees
 maximum safety for the user during handling and avoids any possibility of breakage of the packaging, soft
 packaging and sterilization bags
- The basket is electropolished. This external finish avoids changes in the colour of the basket due to high temperatures during the sterilization cycle and inhibits the formation of water stagnation
- The perforation gap of the base of the basket and the lid is made up of 4x4 mm holes, in order to allow better water penetration during the washing phases and a better drainage capacity. The perforation gap on the sides of the basket is 2x2 mm. The washing process is further facilitated by the fact that the drum is free of shaded areas
- The basket is supplied complete with a lid with a simple and easy closure and a carrying handle. The predisposition for anchoring and fixing the lid is placed on both sides of the basket
- The extremely low weight, thanks to the special production technique, guarantees high comfort during handling and high performance in the drying phase in automatic washing machines
- The design of the silicone instrument holder supports ensures maximum safety and the best handling during the handling of the device
- The simultaneous use of silicone and stainless steel makes it possible to deal with all the processes of decontamination, washing and sterilization with steam and low temperature
- The removable silicone mounts are attached to the body by means of two low-profile stainless steel head screws AISI304
- As for the sterilization process, the basket can be placed in a rigid container or in soft packaging.

Therefore...

- Safe and perfect positioning of microsurgery instruments
- Organize tools quickly and easily
- immediate creation of kits by type of intervention
- use of resistant and high-performance materials (stainless steel and silicone)
- elimination of damage caused to the instrumentation by continuous handling
- Outstanding results after washing








Example with instruments locked with silicone supports ...

EXTREMELY SMALL CONTACT SURFACE FOR EASY CLEANING AND WASHING OF INSTRUMENTS



BEFORE AFTER

SECURE LOCKING OF INSTRUMENTS WITH 2 STAR SILICONE HOLDERS



Legend:

Approximate images and dimensions

The code of the proposed solution does not include the instruments but only the housing system

C.B.M. Srl reserves the right to make any improvements to its products that it deems necessary and useful without compromising the essential characteristics of functionality and safety.











C.B.M. SRL A Socio Unico
Via Castello 10/A
26038 Torre de' Picenardi (Cremona) - Italy
+390375 394 095 - F: +39 0375 394 098
vendite@cbm-srl.com – www.cbm-srl.com

100% Made in Italy

100% Made in CBM